

## MICROKIT STAR PRODUCTS BY SECTORS

### 2.

## INDUSTRIAL WATER, SURFACES AND AIR



(See also bottled and drinking water chapter):



\* KIT P/A *Legionella pneumophila* RPL330, 20 or 200 test or revitalization before use ISO 11731



\* M-Ident-Legionella latex to confirm and serogrouping colonies of *L.pneumophila* KMB301, 50 test



\* Tubes Reductor of Sulfate Desulfobrio American Petroleum Institute Agar TPL048, 20 or 200u



\* RAPIDSCREEN-RAPIDTEST rapid evaluation of microbial bioburden, KAJ001, 20 or 200 test even since 4-8 hours!



\* KITPRO immediate control of hygiene (surfaces and water) by color (changes to blue) KP030, 36u



\* KIT P/A swimming pools 8 microbiological parameters RPL399, 1 or 10 test



\* Dipslides DESINFECTEST® for quantitative control of surfaces and semi-quantitative control in liquids: TPC, yeast & moulds, enterobacteria, coliform, *E.coli*, Staphylococcus, Listeria.



\* MICROFLOW, The most versatile microbiological air sampler UNE 100012 and ISO 14698 : With contact plates, several volumes and several flows.



\* LPT Neutralizing Agar, Total plate count in disinfected surfaces and air, DMT066, 500g, PPL502 Contact Plates UNE 100012



\* Rose Bengal Caf. Agar, Total Yeast & Moulds count avoiding invasion of plates, DMT101, 500g, PPL506 Contact Plates UNE 100012



\* Listeriswabs for searching *L. Monocytogenes* in HACCP MW 570, 125-1000 u